Product Specification				
Balloon Length	: 6, 10, 15 & 20mm			
Balloon Diameter	: 2.25, 2.50, 2.75, 3.00, 3.25, 3.50, 3.75, 4.00,			
	4.50mm			
Balloon Compliance	: Non-Compliant			
Proximal Shaft	: ≤ 0.77mm			
Distal Shaft	: Φ 2.25-4.00mm: ≤ 0.86mm			
	Φ 4.50mm: ≤ 0.93mm			
Usable Length	: 1450mm			
Coating	: Hydrophilic			
Balloon Material	: PET-Nylon 12			
Brachial and Femoral Markers	: 950 and 1050mm			
Guidewire Compatibility	: 0.014"			
Nominal Pressure	: 12 atm			
Rated Burst Pressure	: Φ 2.25-4.00mm: ≥ 22atm, Φ 4.50: ≥ 20atm			



PTCA Balloon Catheter

Cross. Dilate. Optimize

ORDERING INFORMATION

Ø(mm)/Length(mm)	6	10	15	20
2.25	FNC2206	FNC2210	FNC2215	—
2.50	FNC2506	FNC2510	FNC2515	FNC2520
2.75	FNC2706	FNC2710	FNC2715	FNC2720
3.00	FNC3006	FNC3010	FNC3015	FNC3020
3.25	-	FNC3210	FNC3215	FNC3220
3.50	-	FNC3510	FNC3515	FNC3520
3.75	—	FNC3710	FNC3715	FNC3720
4.00	-	FNC4010	FNC4015	FNC4020
4.50	-	FNC4510	FNC4515	FNC4520



Purple Microport Cardiovascular Private Limited 615/616 Avior, L.B.S. Road, Mulund (W), Mumbai-400080

Tel: +(91) 22 25916811/12 Email: info@purplemicroport.com Website : www.purplemicroport.com

All cited trademarks are the property of their respective owners in respective countries. These devices shall only be used for curing diseases on prescription by a licensed physician. Indications contraindications, warnings, potential adverse events and directions for use can be found in instructions for use and product labeling supplied with each device. These devices are prohibited to promote, distribute or sell in any countries or regions where the registration certificates or any other licenses are not obtained. This brochure is not intended for use in aforesaid countries or regions. This brochure is only delivered to specific person for his or her own reference, is not to be distributed widely as an advertisement. Any information contained in this brochure may be changed without prior notice.



