



ORDERING INFORMATION

Product Type	Brand Name	Product Code	Model No	Description	Inflation Device	Y- Connector	Insertion Tool	Torquer	Two/Three way stopcock with Extension line
Kit	Accurex Kit	ACX	ACX	-	√	\checkmark	✓	\checkmark	\checkmark
Inflation Device	X-ACT	ID	AID	X-ACT (20ml, 30atm)	✓				✓
			BID	X-ACT "A" Type (25ml, 30atm)	✓				✓
			CID	X-ACT "S" Type (30ml, 30atm)	✓				✓
			DID	X-ACT "G" Type (20ml, 30atm)	✓				✓
Y-Connector	Y-SYNC	YC	AYC	Y Connector Screw Type		✓			
			BYC	Y Connector Push Pull Type		\checkmark			
			CYC	Y Connector Click Type		✓			
			DYC	Double Y-Connector		✓			
Insertion tool	Insertion Tool	IT	AIT	Insertion tool			✓		
Torquer	Torquer	TR	ATR	Torquer				\checkmark	
Stopcock	Stopcock	SC	ASC	Stopcock_Two way					✓
			BSC	Stopcock_Three way					✓



Purple Microport Cardiovascular Private Limited 615/616 Avior, L.B.S. Road, Mulund (W), Mumbai-400080

Tel: +(91) 22 25916811/12 Email: info@purplemicroport.com Website : www.purplemicroport.com

All cited trademarks are the property of their respective owners in respective countries. These devices shall only be used for curing diseases on prescription by a licensed physician. Indications contraindications, warnings, potential adverse events and directions for use can be found in instructions for use and product labeling supplied with each device. These devices are prohibited to promote, distribute or sell in any countries or regions where the registration certificates or any other licenses are not obtained. This brochure is not intended for use in aforesaid countries or regions. This brochure is only delivered to specific person for his or her own reference, is not to be distributed widely as an advertisement. Any information contained in this brochure may be changed without prior notice.

out prior notice.

ACX/BRO/20210901_03





Inflation Device With Unique Design Features





X-ACT UNIQUE DESIGN FEATURES

- Ergonomic design facilitates excellent eye to hand co-ordination.
- Efficient intra cam locking mechanism allows maintaining pressure and guaranteeing fast and automatic release. Easy deflation even at high pressures.
- Quick latch/push mechanism for gradual and rapid pressure release
- Clear fluid barrel for unobstructed view during bubble purging
- Simple operation and stable performance
- Consistent pressure output
- Single handed release capability for left or right-handed use
- Up to 30atm analog pressure range, provides stable accurate repeatable pressures.
- Available in 20ml (X-Act & X-Act "G" Type), 25ml (X-Act "A" Type) & 30ml (X-Act "S" Type)

INSERTION TOOL:

Guide Wire Insertion Tool is precisely designed for use during PTCA and any other intravascular therapeutic procedures that utilize a guiding catheter.

This provides a mean for inserting guide wire deep into the vasculature.

- Compatible for 0.010" to 0.016" guide wires.
- Needle Length: 100mm
- Inner Diameter: 0.63mm



TORQUER:

Guide Wire Torquer is precisely designed for use PTCA and any other intravascular therapeutic procedures that utilize a guiding catheter.

This provides a mean for maneuvering guide wires once inserted deep into the vasculature.

Compatible for 0.010" to 0.016" guide wires.



Y-CONNECTOR:

Y Connectors are precisely designed for use during PTCA and any other intravascular therapeutic procedures that utilize a guiding catheter.

Guide wire compatibility 7/8/9 F & pressure rating up to 170 psi

Types of Y-connectors:

- a. Y-connector Screw Type
- · Large closing knob makes twisting and valve closing easier
- Minimum blood loss
- b. Y-connector Push Pull Type
- Simple push mechanism allows single handed operation to open & close valve.
- Rapid, efficient & excellent Ergonomics
- · Minimized backflow, limiting blood loss
- c. Y-connector Click Type
- Ergonomic designed provides easy maneuvering and significantly reduced blood loss
- d. Double Y-connector
- · Different valves make it easier to distinguish and operate the two guide wires







