Prod	uct Specification	
Balloon Length	: 8 to 41 mm	
Balloon Diameter	: 1.00, 1.10, 1.25, 1.50 mm	
Balloon Compliance	: Semi Compliant	
Proximal Shaft	: 1.95 F	
Distal Shaft	: 2.70 F	
Usable Length	: 140 cm	
Hydrophilic Coating	: P V P	
Balloon Material	: Nylon Copolymer	
Balloon Folds	: 3 Folds	
Brachial and Femoral Markers	: 90 and 100cm	
Rapid Exchange Length	: 25 cm	
Guide Catheter Compatibility	: 5F	
Guidewire Compatibility	: 0.014"	
Entry Profile	: 0.017"	
Nominal Pressure	: 8 ATM	
Rated Burst Pressure	: 16 ATM, 14 for diameter 3.50, 3.75 with length > 35mm and diameter 4.00 with length > 30mm.	

## **ORDERING INFORMATION**

Length (mm)	1.00	1.10	1.25	1.50
8mm	ACO1.0008	ACO1.1008	ACO1.2508	ACO1.5008
9mm	ACO1.0009	ACO1.1009	ACO1.2509	ACO1.5009
10mm	ACO1.0010	ACO1.1010	ACO1.2510	ACO1.5010
11mm	ACO1.0011	ACO1.1011	AC01.2511	ACO1.5011
12mm	ACO1.0012	ACO1.1012	ACO1.2512	ACO1.5012
13mm	ACO1.0013	ACO1.1013	ACO1.2513	ACO1.5013
14mm	ACO1.0014	ACO1.1014	ACO1.2514	ACO1.5014
15mm	ACO1.0015	ACO1.1015	ACO1.2515	ACO1.5015
16mm	ACO1.0016	ACO1.1016	ACO1.2516	ACO1.5016
17mm	ACO1.0017	ACO1.1017	ACO1.2517	ACO1.5017
18mm	ACO1.0018	ACO1.1018	ACO1.2518	ACO1.5018
19mm	ACO1.0019	ACO1.1019	ACO1.2519	ACO1.5019
20mm	ACO1.0020	ACO1.1020	ACO1.2520	ACO1.5020

Note : Also available in 21 to 41mm lengths

## ⊗ Purple MicroPort

Purple Microport Cardiovascular Private Limited 615/616 Avior, L.B.S. Road, Mulund (W), Mumbai-400080

Tel: +(91) 22 25916811/12 Email: info@purplemicroport.com Website : www.purplemicroport.com

All cited trademarks are the property of their respective owners in respective countries. These devices shall only be used for curing diseases on prescription by a licensed physician. Indications contraindications, warnings, potential adverse events and directions for use can be found in instructions for use and product labeling supplied with each device. These devices are prohibited to promote, distribute or sell in any countries or regions where the registration certificates or any other licenses are not obtained. This brochure is not intended for use in aforesaid countries or regions. This brochure is only delivered to specific person for his or her own reference, is not to be distributed widely as an advertisement. Any information contained in this brochure may be changed without prior notice.

DBCTCO/BRO/20210901\_04





## CTO PTCA Pre-Dilatation Balloon Catheter

