
Product Specification							
Stent Material:	Cobalt Chromium L605						
Polymer:	Biodegradable and Biocompatible polymer (PLGA)						
Drug:	Everolimus						
Drug Dosage:	1.00 µg/mm²						
Strut Thickness:	65 μm (0.065 mm or 0.0026")						
Nominal Pressure:	9 atm						
Rated Burst Pressure:	16 atm, 14 atm with 3.50mm diameter with length > 35 mm and 4.00mm diamete with length >30mm and 4.50 mm diamete						
Mean Recoil:	Average 3.2%						
Mean Foreshortening:	± 1.3 %						







Biodegradable Polymer Based Everolimus Eluting Coronary Stent System

ORDERING INFORMATION

Ø (mm)/ Length (mm)	8mm	13mm	16mm	19mm	24mm	29mm	32mm	37mm	40mm
2.00mm	AC4.1 2.0008	AC4.1 2.0013	AC4.1 2.0016	AC4.1 2.0019	AC4.1 2.0024	AC4.1 2.0029	AC4.1 2.0032	AC4.1 2.0037	AC4.1 2.0040
2.25mm	AC4.1 2.2508	AC4.1 2.2513	AC4.1 2.2516	AC4.1 2.2519	AC4.1 2.2524	AC4.1 2.2529	AC4.1 2.2532	AC4.1 2.2537	AC4.1 2.2540
2.50mm	AC4.1 2.5008	AC4.1 2.5013	AC4.1 2.5016	AC4.1 2.5019	AC4.1 2.5024	AC4.1 2.5029	AC4.1 2.5032	AC4.1 2.5037	AC4.1 2.5040
2.75mm	AC4.1 2.7508	AC4.1 2.7513	AC4.1 2.7516	AC4.1 2.7519	AC4.1 2.7524	AC4.1 2.7529	AC4.1 2.7532	AC4.1 2.7537	AC4.1 2.7540
3.00mm	AC4.1 3.0008	AC4.1 3.0013	AC4.1 3.0016	AC4.1 3.0019	AC4.1 3.0024	AC4.1 3.0029	AC4.1 3.0032	AC4.1 3.0037	AC4.1 3.0040
3.5mm	AC4.1 3.5008	AC4.1 3.5013	AC4.1 3.5016	AC4.1 3.5019	AC4.1 3.5024	AC4.1 3.5029	AC4.1 3.5032	AC4.1 3.5037	AC4.1 3.5040
4.00mm	AC4.1 4.0008	AC4.1 4.0013	AC4.1 4.0016	AC4.1 4.0019	AC4.1 4.0024	AC4.1 4.0029	AC4.1 4.0032	AC41 4.0037	AC4.1 4.0040
4.50mm	AC4.1 4.5008	AC4.1 4.5013	AC4.1 4.5016	AC4.1 4.5019	AC4.1 4.5024	AC4.1 4.5029	AC4.1 4.5032	AC4.1 4.5037	AC4.1 4.5040

Note : Also available in 44 & 48mm lengths



Purple Microport Cardiovascular Private Limited 615/616 Avior, L.B.S. Road, Mulund (W), Mumbai-400080

Tel: +(91) 22 25916811/12 Email: info@purplemicroport.com Website : www.purplemicroport.com

All cited trademarks are the property of their respective owners in respective countries. These devices shall only be used for curing diseases on prescription by a licensed physician. Indications contraindications, warnings, potential adverse events and directions for use can be found in instructions for use and product labeling supplied with each device. These devices are prohibited to promote, distribute or sell in any countries or regions where the registration certificates or any other licenses are not obtained. This brochure is not intended for use in aforesaid countries or regions. This brochure is only delivered to specific person for his or her own reference, is not to be distributed widely as an advertisement. Any information contained in this brochure may be changed without prior notice.

